


SAC2 Adhesive Conductive Paste

SAC2 is an adhesive EEG electrode cream to use with surface electrodes. SAC2 is intended to be applied on intact skin. SAC2 is intended for use in clinical and research EEG/EP recordings from humans.

Item code:	Item Description:	Picture:
SAC2-10	SAC2 Adhesive Conductive Paste	

Delivery time: Within 5 days

Term of Payment: 30 Days after order placed

It is used with external electrodes as the conductor between the scalp and recessed electrodes to reduce impedance between the electrode surface and the skin.

Features:

- High adhesion
- Viscous and white cream
- PH: 8.5

Packaging

Quantity of each box	10 Tubes (100g per tube)
Sterilization	NO
Operation ambient temperature	15°C - 30°C
Relative humidity	Protect against sunlight
Storage temperature	15°C - 30°C