



Description

The Dragonfly laryngeal electrodes are made by flexible material to avoid any risk of injury during the positioning of the endotracheal tube. Dragonfly can be applied on any endotracheal tube armed and not armed.

UMDNS 11441

Materials

Safety

Specifications

Environment

Electrode	Ag
Substrate	PET
Lead wire	Cu/Sn with PVC jacket
Connector	CuZn with PVC jacket
Packaging	Tyvek/PET-PE

Biocompatibility (UNI EN ISO 10993)	Yes
Sterilization method	Ethylene oxide (EtO)
Latex	No
Phthalates	No

Generic codes	LSEXXXXXXXXXX
Number of channels	1 - 2
Compatible ET tube diameter	2.0-5.5 mm (pediatric version) 6.0-10.0 mm (adult version)
Cable length	75 - 200 cm
Connector	Touch Proof 1,5 mm DIN 42802
Cable colour	● ●

RoHS II	✓
REACH	✓